


Analyzing Fusion Genes with Next-Generation Sequencing Technology




Raed Samara, Ph.D
Global Product Manager



- QIAGEN products shown here are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.
- For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at **www.qiagen.com** or can be requested from QIAGEN Technical Services or your local distributor.

- 1 Fusion genes: What they are and their historical perspective
- 2 Fusion gene detection: Current status
- 3 RNA sequencing vs. digital RNA sequencing
- 4 How to detect and accurately quantify novel fusion genes in your sample


- 1 Fusion genes: What they are and their historical perspective
- 2 Fusion gene detection: Current status
- 3 RNA sequencing vs. digital RNA sequencing
- 4 How to detect and accurately quantify novel fusion genes in your sample



What is a gene fusion?

A new gene formed by the chromosomal translocation of two parts of different genes

A. Chromosomal Translocation




Source: Leonard, G. [CC BY-SA 3.0 (<http://creativecommons.org/licenses/by-sa/3.0/>)], via Wikimedia Commons
https://upload.wikimedia.org/wikipedia/commons/b/b0/Gene_Fusion_Types.png

Why are gene fusions important?

They are implicated in hematologic and soft tissue malignancies...

Timeline of gene fusion discoveries in:
Hematologic malignancies
Soft tissue malignancies



Source: Kumar-Sinha, C., Kalyana-Sundaram, S., Chinnaiyan, A.M. (2015) Landscape of gene fusions in epithelial cancers: seq and ye shall find. Genome Med. 7, 129.


Why are gene fusions important?

And in rare and common epithelial cancer...

Timeline of gene fusion discoveries in:


Rare epithelial cancers

Common epithelial cancers




Source: Kumar-Sinha, C., Kalyana-Sundaram, S., Chinnaiyan, A.M. (2015) Landscape of gene fusions in epithelial cancers: seq and ye shall find. Genome Med. 7, 129.

- 1 Fusion genes: What they are and their historical perspective
- 2 Fusion gene detection: Current status**
- 3 RNA sequencing vs. digital RNA sequencing
- 4 How to detect and accurately quantify novel fusion genes in your sample



At an unprecedented rapid pace...




Source: Mertens, F., Johansson, B., Fioretos, T., Mitelman, F. (2015) The emerging complexity of gene fusions in cancer. Nat Rev Cancer 15, 371–81.

Current technologies: Advantages and disadvantages

Current technology	Advantages	Disadvantages
PCR-based	Accuracy	Limited sample and assay throughput Requires a lot of RNA
Whole transcriptome sequencing (WTS)	Throughput power	Expensive Complex data
FISH (Fluorescence in situ hybridization)	Established method for routine testing	Laborious
Conventional Targeted RNA sequencing	Manageable data Relatively low per-sample cost	PCR duplicates, which limit the sensitivity of the panel and reduce the ability to confidently detect low-abundance fusions

- 1 Fusion genes: What they are and their historical perspective
- 2 Fusion gene detection: Current status
- 3 RNA sequencing vs. digital RNA sequencing**
- 4 How to detect and accurately quantify novel fusion genes in your sample

PCR amplification is required for target enrichment, but...



PCR duplicates limit accurate quantification

Conventional
targeted RNA
sequencing



EGFR exon 21

- Five library fragments that look exactly the same.
Cannot tell whether they represent:
1. Five **unique** fusion transcripts, or
 2. Five **copies** of the **same** fusion transcript (**PCR duplicates**)




Quantification based on non-unique reads does
not reflect quantities of original fusion molecules



If you cannot eliminate PCR duplicates, correct for them

Capture and freeze original status of fusion transcripts by tagging them before any amplification



Tag (barcode) to identify unique RNA molecules

TATCGTACAGAT
(12 nucleotides long)

Incorporate this **random** barcode (signature) into the original RNA molecules (**before** amplification) to preserve their uniqueness

Count and analyze single original molecules (not total reads) = digital sequencing

Conventional
targeted RNA
sequencing



Five library fragments that look exactly the same.
Cannot tell whether they represent:


1. Five **unique** fusion transcripts, or
2. Five **copies** of the **same** fusion transcript (**PCR duplicates**)

UMIs **before** any
amplification

Digital sequencing
with UMIs




Five **unique** fusion transcripts
since 5 molecular barcodes are detected




Five **copies** of the **same** fusion transcript (**PCR duplicates**)
since 1 molecular barcode is detected

- 1 Fusion genes: What they are and their historical perspective
- 2 Fusion gene detection: Current status
- 3 RNA sequencing vs. digital RNA sequencing
- 4 How to detect and accurately quantify novel fusion genes in your sample**

QIAseq targeted RNAscan workflow




Panels, molecularly-barcoded (UMI) adapters and data analysis algorithms



QIAGEN's solution: QIAseq targeted RNAscan Panels

Current technology	Disadvantages	How QIAseq targeted RNAscan panels address disadvantages of current technologies
PCR-based	Limited sample and assay throughput Requires a lot of RNA	Profile hundreds of gene fusions in up to 384 samples simultaneously Requires 10 ng RNA
Whole transcriptome sequencing (WTS)	Expensive Complex data	Cost-effective Easy data analysis
FISH (Fluorescence in situ hybridization)	Laborious	Streamlined, automation-friendly protocol
Conventional Targeted RNA sequencing	PCR duplicates, which reduce the ability to confidently detect low-abundance fusions	Digital sequencing removes PCR duplicates to increase sensitivity of panel, and detect low-abundance fusions


As low as 10 ng



RNA conversion step is crucial


More RNA converted \longrightarrow More fusions detected

- Test kits:
 - QIAseq Targeted RNAscan panel kit
 - **Competitor A** panel kit
- Sample: For each reaction, 10 ng universal RNA spiked with 1 pg of an artificial RNA molecule.
- Follow each kit's standard workflow to generate 2nd strand product.
- Detect converted spike-in RNA with qPCR
- Result showed RNAscan has superior RNA conversion efficiency.
- Higher conversion efficiency means higher sensitivity, especially with low input.




	RNAscan	A
C_T (average \pm SD)	22.8 \pm 0.05	38.5 \pm 0.53

As low as 10 ng



QIAseq Targeted RNAscan Panel: Workflow

As low as 10 ng



Count and analyze single original molecules (not total reads) = digital sequencing

Conventional targeted RNA sequencing




Five library fragments that look exactly the same.
 Cannot tell whether they represent:
 1. Five **unique** fusion transcripts, or
 2. Five **copies** of the **same** fusion transcript (PCR duplicates)

Molecular barcodes
 before any amplification

Digital sequencing with UMIs




Five **unique** fusion transcripts since 5 molecular barcodes are detected



Five **copies** of the **same** fusion transcript (PCR duplicates) since 1 molecular barcode is detected


As low as 10 ng



Up to 1000 primers
in a single pool

QIAseq Targeted RNAscan Panel: Workflow (L)

As low as 10 ng



Up to 1000 primers
in a single pool


Only one primer is
needed

SPE



No need for long fragments to give space for two primers apart from the fusion junction, so shorter fragments are better tolerated. This approach will help deal with FFPE sample challenges where many RNA fragments may not have enough length to cover both GS1 and GS2 sequence in competitor A panels. Competitor A panels might require longer fragments, which will reduce the efficiency of PCR and the crossing of spanning reads across fusion junctions.

Nested PCR



One gene/target-specific primer

Two gene/target-specific primers



- Flexibility in primer design
- Generate relatively small library fragments to maintain compatibility with fragmented RNA (FFPE samples)

Design based on junction


Primers target both ends of a fusion




fusion pair	5' chr	start	end	strand	3' chr	start	end	strand	ENSTs and exons
ETV6-NTRK3	12	11853560	11853561	+	15	87940752	87940753	-	ETV6 (ENST00000396373) exon 4 to NTRK3 (ENST00000394480) exon 15

Exon- or gene-based design for unknown partners

Primers target one end of a fusion



As low as 10 ng




- Up to 1000 primers in a single pool
- Dual sample indexing system; up to 384 samples

QIAseq Targeted RNAscan Panel: Workflow (I)




8–9 hours



QIAseq Targeted RNAscan Panel: Workflow (L)



8–9 hours



Feature	Benefit
Low RNA input (as low as 10 ng of unenriched RNA)	Preserve sample
Unique RNA conversion chemistry	Convert more RNA to cDNA to detect more fusions
Unique molecular indices (UMIs)	Enables single molecule counting, removal of PCR duplicates and error correction; these features correct for amplification bias, distinguish clonal replicates and deliver absolute measurements of fusion transcripts
High primer multiplexing capability (up to 1000 primers)	Detects a large number of fusions
Single pool of primers	Easier sample handling
SPE-based target enrichment	Flexibility in primer design Generate relatively small library fragments to maintain compatibility with fragmented RNA (FFPE samples)
High sample multiplexing (up to 384 samples)	Increased sample throughput to decrease sequencing costs
No prior knowledge of breakpoint needed	Discover novel partners
Automation-friendly workflow	Streamline operations for high throughput
Suite of complementary data analysis tools	Save resources
Very cost-effective compared to competition	Save \$\$\$

Two boxes are needed to support the workflow

Panel box (kit)



- Reagents for cDNA synthesis
- Fusion-specific primers for enrichment
- All required buffers and enzymes
- Magnetic beads

Index box (kit)



- Molecularly-barcoded library adapters, primers to prepare libraries that are sample indexed and sequencing platform-specific

At a glance

Panels

- RNAscan
 - Cataloged
 - Custom
 - Extended

Indices

- Illumina
 - 12-index
 - 96-index (4 sets, for up to 384-plex)
- Ion Torrent
 - 12-index
 - 96-index

List of panels

Panel	Variant (cat) number	Number of primers	Types of coverage
Hematology panel	FHS-001Z	156	Breakpoint
Solid tumor panel	FHS-002Z	101	Breakpoint
Lung cancer panel	FHS-003Z	137	Breakpoint
Oncology panel	FHS-3001Z	950	Breakpoint

Primers target both ends of a fusion



For QIAseq targeted RNAscan panels

Cat ID	Name	Sequencing platform	Number of samples per sequencing run that can be multiplexed	Number of samples each kit (SAP ID) can process
333714	QIAseq 12-Index I (48)	Illumina	12 samples per sequencing run	48 samples
333727	QIAseq 96-Index I Set A (384)	Illumina	96 samples per sequencing run	384 samples
333737	QIAseq 96-Index I Set B (384)	Illumina	96 samples per sequencing run (192 if used with Set A)	384 samples
333747	QIAseq 96-Index I Set C (384)	Illumina	96 samples per sequencing run (288 if used with Sets A, B)	384 samples
333757	QIAseq 96-Index I Set D (384)	Illumina	96 samples per sequencing run (384 if used with Sets A, B, C)	384 samples
333764	QIAseq 12-Index L (48)	Ion Torrent	12 samples per sequencing run	48 samples
333777	QIAseq 96-Index L (384)	Ion Torrent	96 samples per sequencing run	384 samples

Application: Proof of sensitivity and specificity of the QIAseq Targeted RNAscan panels

- Engineered positive control for 12 popular fusions as testing material
- QIAseq Targeted RNAscan panel was designed based on fusion's breakpoint and sequence used from both genes
- Control FFPE sample was used with 15 ng as starting material; same amount of normal universal RNA was used as negative control
- Library was run with sequencer MiSeq
- Fusion was detected with current QIAseq Targeted RNAscan data analysis pipeline

Fusion Partners included

	RNA Fusion	Primary Cancer Tissue
1	EML4-ALK	Lung
2	NPM1-ALK	Lymphoid
3	KIF5B-RET	Lung
4	NCOA4-RET	Thyroid
5	CD74-ROS1	Lung
6	SLC34A2-ROS1	Lung, Stomach
7	TPM3-NTRK1	Lung, Large Intestine
8	TFG-NTRK1	Thyroid (rare)
9	FGFR3-BAIAP2L1	Urinary tract (rare)
10	FGFR3-TACC3	Urinary tract, CNS
11	PAX8-PPARG	Thyroid
12	ETV6-NTRK3	Kidney, Breast, Soft Tissue

100% sensitivity (12/12 positive), 0% positive

Fusion in sample	UMIs for detected crossing-junction reads
EML4-ALK	18
NPM1-ALK	5
KIF5B-RET	25
NCOA4-RET	46
CD74-ROS1	11
SLC34A2-ROS1	8
TPM3-NTRK1	46
TFG-NTRK1	5
FGFR3-BAIAP2L1	37
FGFR3-TACC3	60
PAX8-PPARG	28
ETV6-NTRK3	30

● All fusions in control were detected with at least 5 MTs

Example of positive fusion detected: In positive controls, FGFR3-TACC3 fusion is one the 12 positive fusions

FGFR3--TACC3

Summary

Role	Gene(s)	Donor			Acceptor			Links	Crossing		Spanning	
		Chr	Position	Strand	Chr	Position	Strand		Reads	Tags	Reads	Tags
Donor	FGFR3	chr4	1,806,935	+	chr4	1,807,115	+	1	32	15	0	0
Acceptor	TACC3	chr4	1,739,701	+	chr4	1,739,701	+	2	720	459	0	0
Fusion	FGFR3--TACC3	chr4	1,806,935	+	chr4	1,739,701	+	1	99	60	0	0

Confidently
detect low-
abundance
fusions

Molecular barcodes

- Correct for PCR duplicates
- Unmatched sensitivity

Flexibility to
detect known
and novel
fusions

Single primer extension and powerful pipelines

- Only one of the two targets need to be known/targeted
- Targets are defined based on breakpoint, exons or genes

Completely integrated and optimized workflow: From Sample to Insight

Digital sequencing with molecular barcodes

- Remove PCR duplicates and errors for accurate quantification

Barcode-aware, comprehensive data analysis pipelines

- Confidently call low-abundance fusion transcripts and identify novel variants

Single primer extension approach

- No need for nested PCR primer design; ability to detect more fusions

Throughput flexibility

- Molecularly profile hundreds of fusion transcripts in hundreds of samples, simultaneously

Compatibility with fragmented or preserved RNA

- Sequence RNA from FFPE samples or liquid biopsy samples

Fast custom design and manufacturing turnaround time

- Receive custom panel or extend an existing panel in 2–3 weeks.

- Barcode-aware fusion caller has been developed
- Caller is available on the cloud
- In conjunction with molecular barcodes incorporated in the workflow, the caller can confidently call low-abundance fusions (down to 1% expression level)
- Fusion caller will do the following:
 - Mapping
 - Alignment
 - Molecular barcode counting
 - Fusion calling
 - Fusion annotation – based on internal curated database

Inputs

FASTQ or BAM files are uploaded into cloud-based data analysis portal


The following inputs are needed (by user):

- Set up file
- Panel used
- File lanes
 - 1-lane (MiSeq/HiSeq/NextSeq concatenated)
 - 4-lane (NextSeq individual lane files)


Outputs


Summary file

- Stats
 - Specificity
 - Uniformity
 - Molecular barcode counts
- Fusions
 - Curated/novel/no reference available
 - Statistical confidence
 - Graphical representations



Different panels can be plugged into the same targeted NGS workflow







Questions?



Contact QIAGEN

Call: 1-800-426-8157

Email:

techservice-na@QIAGEN.com

BRCsupport@QIAGEN.com

QIAseq.NGS@QIAGEN.com

QIAWebinars@QIAGEN.com