



## Targeted RNAseq for gene expression using unique molecular indexes (UMIs): Introduction to QIAseq Targeted RNA Panels.

Samuel Rulli, Ph. D, Global Product Manager QIAseq Targeted RNA Panels



- QIAGEN products shown here are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.
- For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.QIAGEN.com](http://www.QIAGEN.com) or can be requested from QIAGEN Technical Services or your local distributor.



- 1 Using NGS approaches for gene expression analysis
- 2 Principle of QIASEQ Targeted RNAseq
  - 2.1 Molecular Barcodes
  - 2.2 QIASEQ RNA workflow
- 3 An application of the QIASEQ RNA system
  - 3.1 QIASEQ data analysis
  - 3.2 Ingenuity IPA
- 4 Uses of UMIs in other types of targeted sequencing
- 5 Summary and Discussion




- 1 Using NGS approaches for gene expression analysis
- 2 Principle of QIAseq Targeted RNAseq
  - 2.1 Molecular Barcodes
  - 2.2 QIAseq RNA workflow
- 3 An application of the QIAseq RNA system
  - 3.1 QIAseq data analysis
  - 3.2 Ingenuity IPA
- 4 Uses of UMIs in other types of targeted sequencing
- 5 Summary and Discussion

# Are “reads” the best we can do with sequencing?

## Whole transcriptome sequencing for gene expression

- Quantifies and characterizes all RNA
- Final data point in READS per target






## Whole transcriptome sequencing for gene expression

- Quantifies and characterizes all RNA
- Final data point in READS per target

## Drawbacks

- Complex library construction
- 1µg of total RNA
- Fails on FFPE & fragmented RNA
- Large computational requirements
  - Massive amount of data generated
  - Filtering, alignment, assembly, curation
  - Aggressive normalization for quantification
  - Not at all straightforward
  - Requires skilled bioinformatics scientists
- Cost
  - Large read budget = money
  - Limits sample numbers in studies
- Only runs on HT instruments
  - Limits accessibility to core labs



What are the advantages of applying **targeted gene profiling** to NGS?

- Use read budget only for genes of interest
  - Cost
  - Time (quick prep, run, analysis)
  - Sample throughput - multiplex many samples
- Desktop platforms can now be used for RNA analysis (500,000 reads per sample instead of 20,000,000 or more)
- Simplified bioinformatics (no assembly required)
  - Don't need that bioinformatics guy down the hall
- Minimal sample pre-processing
  - No ribosomal depletion or blocking
  - No polyA selection
  - Only nanogram quantities of RNA required
  - 6 hour sample prep - only need thermocycler and magnet.

When? Who?

- Scientists with known gene list or pathway
- Follow up on WTS or microarray
- Alternative to digital PCR/ Nanostring/Taqman qPCR/ Fluidigm/OpenArray/Wafergen qPCR


Targeted RNAseq still is a “read” based approach to understanding gene expression.

How do we go from “reads” to counting transcripts?



PCR duplicates and amplification bias are major issues in current RNAseq workflows, as they result in biased and inaccurate gene expression profiles

- PCR duplicates
- Amplification bias
- Ratio based on reads



Molecular barcodes each “capture” event

12 random bases  
16 million barcodes

5' AATGTACAGTATTGCGTTTTGNNNNNNNNNNNNNNCGGCAGGAGACGAAGAG 3'

RS2

MT

GS

Use approximately 1/20<sup>th</sup> of total number of primers in a reaction - each gene specific primer is statistically unique!





- 1 Using NGS approaches for gene expression analysis
- 2 Principle of QIASEQ Targeted RNAseq
  - 2.1 Molecular Barcodes
  - 2.2 QIASEQ RNA workflow
- 3 An application of the QIASEQ RNA system
  - 3.1 QIASEQ data analysis
  - 3.2 Ingenuity IPA
- 4 Uses of UMIs in other types of targeted sequencing
- 5 Summary and Discussion

Molecular barcodes allow the counting of original gene levels instead of PCR duplicates, thereby enabling digital sequencing and resulting in unbiased and accurate gene expression profiles


Tag each gene with unique molecular barcodes

Count unique barcodes, not reads







“Leave no scientist behind” ....








Everything needed to go from RNA → Library in one kit, one day!



# 'Simple' procedure, integrated library preparation



## Customer Criteria

## Differential Gene Expression By QIAseq NGS

Species coverage

Human - Catalog, Extended, Virtual, Custom panels  
Mouse + Rat - Custom

Biological replicates

Essential for robustness of experimental design (and statistics!)

Short reads for FFPE, and Exosomal RNA

Average amplicon 97 bps'; range 95-130 bases

Coverage across the transcript (i.e. cover every exon)

We are counting single common regions per gene. Same design philosophy as RT2 PCR Arrays

Depth of sequencing

High enough to infer accurate statistics  
*determined by Smcounter - ~2-5 reads per random barcode*

Role of sequencing depth

Capture enough unique tags of each transcript such that statistical inferences can be made (>10 tags per gene)

Stranded library prep



Not required, amplicons do not overlap lncRNA

Type of reads (paired or Unpaired?)

Not necessary, 150 base single reads more than enough for accurate data


mRNA and lncRNAs

QIAseq was designed against database containing lncRNA and mRNA. Assay are specific for lncRNA or mRNA. Currently 54,881 genes from Ensembl version 81



## Free circulating nucleic acids

RNA and DNA from dead cells shed into the bloodstream, can contain cancer-related mutations.




## Exosomes

Tiny microvesicles found in body fluids that transport RNA between cells.



## Circulating tumor cells


Tumor cells shed from a tumor into the bloodstream carrying genetic information.



## Tissue samples

Fresh tissue or archived **FFPE** samples

● QIAGEN comprehensive sample isolation portfolio compatible with QIAseq RNA



A liquid biopsy is a liquid biomarker that can be isolated from body fluids, such as blood, saliva, urine, ascites, or pleural effusion. Like a tissue biopsy, it is a representative of the tissue from which it has spread.


Liquid biopsies have become more clinically useful in recent years due to the ability to pair tests on circulating tumor cells with genomic tests.

Diaz, Jr., L.A. and Bardelli, A. (2014) "Liquid biopsies: genotyping circulating tumor DNA." *Am. Soc. Clin. Oncol.* **32**, 579.




# Liquid biopsy: Circulating biomarkers for cancer

Tumors shed both intact cells (resulting in circulating tumor cells) as well as cellular components, such as nucleic acids (resulting in cell-free DNA or RNA).



Exosomes: Small membrane vesicles (30–100 nm), secreted by most cell types into the bloodstream.



## Functional biomolecules:

- DNA fragments (exosomal DNA, exoDNA)
- Proteins and/or peptides
- mRNA
- microRNA (miRNA)
- Lipids


- Exosomes play a central role in cell-to-cell communication
- The majority of DNA associated with tumor exosomes is double-stranded, representing whole genomic DNA
- Biological molecules (protein, RNA, and miRNA) contained in exosomes are well protected by a lipid bilayer membrane that confers a high degree of stability

Rolfo, C. et al. (2014) "Liquid biopsies in lung cancer: the new ambrosia of researchers." *Biochimica et Biophysica Acta* **1846**, 539.

Klevebring, D. et. al. (2014) "Evaluation of exome sequencing to estimate tumor burden in plasma." *PLOS One* **9**, e104417.

From sample to extracellular vesicle RNA isolation in just 1 hour


## Separate Serum/Plasma



## Isolate Exosomes




## Isolate RNA




- Microvesicle isolation
  - 20 minutes
- RNA isolation
  - 35 minutes

Co-precipitation of large protein complexes

## Ultracentrifugation (UC)




## Eluate from exoEasy



- Both preparations contain vesicle-shaped structures within an expected size range
- **UC:** Many smaller, unidentified structures/particles that do not match the expected size
- **exoEasy:** Intact vesicles with higher purity

## Bioanalyzer sizing



2 ml plasma was pre-filtered (0.8  $\mu\text{m}$ ) to exclude larger particles.

Both methods purify RNA of similar size and yield



# EVs contain both small and large RNAs

Purification of large, intact, non-degraded RNAs from EVs

exoRNeasy total RNA


Small RNA fraction

Large RNA fraction






mRNA exclusively within vesicles – near 100% bound



miRNA in vesicular and non-vesicular fractions  
(e.g. free Ago2 complexes)


\* Arroyo, J.D. et al. (2011) Argonaute2 complexes carry a population of circulating microRNAs independent of vesicles in human plasma. *Proc. Natl. Acad. Sci. USA* **108**, 5003.



Flexible experiment design for any Catalog Panel options:



## Comprehensive Panels (available for 12, 96 or 384 samples)

- Cancer Transcriptome (395)
- Inflammation & Immunity Transcriptome (475)
- Signal Transduction PathwayFinder (406)
- Stem Cell & Differentiation Markers (293)
- Molecular Toxicology Transcriptome (370)
- Angiogenesis & Endothelial Cell Biology (340)
- Apoptosis & Cell Death (264)
- ECM & Adhesion Molecules (421)



Flexible experiment design for any Extended Panel options:

Add 25 of your favorite targets (mRNAs or lncRNAs) to QIAGEN's comprehensive panel




What is the role of tumor suppressor lncRNAs?  
Find out!




Flexible experiment design for any pathway panel options:

## QIaseq Targeted RNA pathway panels (available for 12, 96 or 384 samples)


Each panel contains 84 genes + controls and housekeeping genes.  
Choose from over 180 panels!



Pathways




Diseases




miRNA Targets

## Built-in controls:




- gDNA assays to control for any gDNA contamination in the RNA sample
- Mean tags per target calculated and mRNAs near this number are flagged during analysis as 'close to noise level'



- Multiple HKG assays to normalize data to make sample-to-sample and run-to-run comparisons possible
- Flexible – use one, two, all, none or any other genes as normalizers
- HKG efficacy evaluation built into secondary data analysis

## Online custom builder:



Design for gene | View design requests

Enter or upload your search list: ?


Enter search terms

Help and search term examples

Upload a file

Human (Homo sapiens)

Search



Panel Details

Add controls

GDC controls (6 as a set)

select house keeping genes (HKGs)

all house keeping genes (HKGs) - (10 as a set)

next


- Choose your own gene content from 54,881 human genes and lncRNA
- Easy to use online Custom Panel Builder to tailor panel to your research needs
  - Input list of genes
  - Select proper controls (genomic DNA contamination control, HKGs, or your own)
  - Output: list of genomic coordinates for primers designed specifically for your genes of interest



Custom panel number

All your custom designs are saved for easy retrieval

Have questions? Easily contact us



The screenshot shows a web interface for managing custom designs. At the top right, there is a 'Contact us for any questions' link. Below this is a navigation bar with 'Design for gene' and 'View design requests'. A table lists design requests with columns for 'Request ID' and 'Details'. One request is visible with ID '16009', 'GeneGlobe cat#: CRHS-10034Z-29', and 'Date designed: February 19, 2016 - 11:53 PM CET'. A 'Download' link is present in the details. A 'configure' button is located at the bottom right of the table area.

Request ID	Details
16009	GeneGlobe cat#: CRHS-10034Z-29 Date designed: February 19, 2016 - 11:53 PM CET Download

Download zip file containing:

- Summary file
- Bed file

Configure and order

Gene ID and symbol

Strand of the genome the gene is on

Amplicon coordinates


Designated controls are shown here

	A	B	C	D	E	F	G	H	I	J	K	L	M
1	gene id	gene symbol	gene strand	chrom	GRCh38 loc 5	GRCh38 loc 3	control type	single exon	# Gencode Basic RNAs	# Gencode Basic RNAs matched	# off target genes	off target genes	amplicon not genome unique
2	ENSG00000041988	THAP3	1	chr1	6628529	6628647	reference_gene	1	3	3	0		0
3	ENSG00000084072	PPIE	1	chr1	39743264	39743891	reference_gene	0	4	4	0		0
4	_GDC_CONTROL_07_	_GDC_CONTROL_07_	1	chr1	104793033	104793127	gDNA_control	1	0	0	0		0
5	ENSG00000213281	NRAS	-1	chr1	114713799	114713893		1	1	1	0		0
6	ENSG00000174775	HRAS	-1	chr11	532251	532355		1	5	4	0		0

- Single exon (1) means both primers are within one exon
- **# Gencode basic RNAs:** Total number of RNA transcripts found for the gene in Gencode
- **# Gencode basic RNAs matched:** # of RNA transcripts targeted by the designed amplicon
- **# off target genes:** Rough prediction of # of off target genes that will also get enriched by the primer pair for the target gene
- **Amplicon not genome unique:** Reads that will not be able to be uniquely mapped to the genome; Some MT counts might come from another loci

## Bed file


Location of designed amplicon






- 1 Using NGS approaches for gene expression analysis
- 2 Principle of QIASEQ Targeted RNAseq
  - 2.1 Molecular Barcodes
  - 2.2 QIASEQ RNA workflow
  - 2.3 QIASEQ RNA performance
- 3 An application of the QIASEQ RNA system**
  - 3.1 QIASEQ data analysis
  - 3.2 Ingenuity IPA
- 4 Uses of UMIs in other types of targeted sequencing
- 5 Summary and Discussion

## Small Molecules – Signal Transduction Application




- HEK293T Cells were treated with 90 different chemical inhibitors
- The 421 Signal Transduction Gene QIAseq Panel was interrogated
- In one day we went from total RNA to sequence ready libraries for 96 samples. The final libraries were quantified, normalized, and pooled. Prior to loading onto a NextSeq, the denatured libraries were diluted to the appropriate input concentration to obtain to generate suitable clusters on the NextSeq



Cells in  
96 well  
plates



## Purification



RNA isolation  
of 96  
samples

## Quality control

RNA Integrity  
(96 samples  
done  
automatically  
while at  
lunch)





RNA  
quantification  
(16 samples/  
90 secs)



## Small Molecules – Signal Transduction Application


- HEK293T Cells were treated with 90 different chemical inhibitors
- The 421 Signal Transduction Gene QIAseq Panel was interrogated
- In one day we went from total RNA to sequence ready libraries for 96 samples. The final libraries were quantified, normalized, and pooled. Prior to loading onto a NextSeq, the denatured libraries were diluted to the appropriate input concentration to obtain to generate suitable clusters on the NextSeq





## Small Molecules – Signal Transduction Application


- HEK293T Cells were treated with 90 different chemical inhibitors
- The 421 Signal Transduction Gene QIAseq Panel was interrogated
- In one day we went from total RNA to sequence ready libraries for 96 samples. The final libraries were quantified, normalized, and pooled. Prior to loading onto a NextSeq, the denatured libraries were diluted to the appropriate input concentration to obtain to generate suitable clusters on the NextSeq
- The parameters of the NextSeq sequencing run were; Single 151 bp read, with a Custom Sequencing Primer (included in kit)



Normalized, pooled libraries



## QIAseq Targeted RNA Data Analysis automated workflow



- **Read Mapping**

- Identify the possible position of the read within the reference genome
- Align the read sequence to reference sequences

- **Primer Trimming**

- Remove the primer sequences from the reads

- **Molecular Barcode Counting**

Go get coffee



## Primary data analysis for QIAseq Targeted RNA Sequencing

	A	B	C	D	E	F	G	H	I	J	K	L
1	read set	16693.ADCY_5_S67	16693.AKT_S89	16693.AMP_K_S3	16693.ATM_S36	16693.AURO_RA_S4	16693.BCR-ABL_S41	16693.Ca2-ATPase_S73	16693.Ca2-channel_S77	16693.CaMK-II_S6	16693.CASP_1_S50	16693.CASP3-6-7-8-10_S52
2	reads total	2,121,847	903,510	2,339,304	1,963,219	2,213,291	1,309,441	953,046	2,132,040	2,303,187	1,212,167	1,948,387
3	reads dropped, < 55 bp	38,535	46,417	79,134	55,341	122,499	16,812	57,348	48,492	40,956	14,611	19,011
4	reads dropped, MT has >1 N base	0	0	0	0	0	0	0	0	0	0	0
5	reads dropped, not mapped	6,577	4,676	13,526	8,591	14,814	4,149	5,595	7,834	6,528	2,940	5,411
6	reads dropped, < 60 bp aligned	51	31	137	48	163	40	47	80	65	27	71
7	reads dropped, off-target	2,029	1,623	6,819	2,726	5,097	1,482	1,597	2,313	2,033	1,185	2,141
8	reads dropped, poor alignment near MT	62,481	26,200	69,538	62,415	65,313	42,108	29,258	64,636	70,070	37,000	59,771
9	reads dropped, high reads per MT	969,174	399,563	1,664,150	1,199,098	1,679,405	0	562,201	1,045,685	637,535	0	0
10	reads used for MT counting	1,043,000	425,000	506,000	635,000	326,000	1,244,850	297,000	963,000	1,546,000	1,156,404	1,861,967
11	MTs	68,726	27,897	35,963	41,171	22,704	83,902	20,957	63,003	103,094	106,252	126,204
12	reads per MT mean	15.2	15.2	14.1	15.4	14.4	14.8	14.2	15.3	15.0	10.9	14.1
13												
14	reads dropped, < 55 bp	1.8%	5.1%	3.4%	2.8%	5.5%	1.3%	6.0%	2.3%	1.8%	1.2%	1.0%
15	reads dropped, MT has >1 N base	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
16	reads dropped, not mapped	0.3%	0.5%	0.6%	0.4%	0.7%	0.3%	0.6%	0.4%	0.3%	0.2%	0.3%
17	reads dropped, < 60 bp aligned	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
18	reads dropped, off-target	0.1%	0.2%	0.3%	0.1%	0.2%	0.1%	0.2%	0.1%	0.1%	0.1%	0.1%
19	reads dropped, poor alignment near MT	2.9%	2.9%	3.0%	3.2%	3.0%	3.2%	3.1%	3.0%	3.0%	3.1%	3.1%
20	reads dropped, high reads per MT	45.7%	44.2%	71.1%	61.1%	75.9%	0.0%	59.0%	49.1%	27.7%	0.0%	0.0%
21	reads used for MT counting	49.2%	47.0%	21.6%	32.3%	14.7%	95.1%	31.2%	45.2%	67.1%	95.4%	95.6%

## Primary data analysis for QIAseq Targeted RNA Sequencing

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T
1	gene id	gene symbol	gene strand	chrom	loc 5'	loc 3'	control type	single exon	16693.ADCY_5_S67	16693.AKT_S89	16693.AMP_K_S3	16693.ATM_S36	16693.AURO_RA_S4	16693.BCR-ABL_S41	16693.Ca2-ATPase_S73	16693.Ca2-channel_S7	16693.CaMK-II_S6	16693.CASP_1_S50	16693.CASP_3-6-7-8-10_S52	16693.CDK_S9
397	ENSG00000115596	WNT6	1	chr2	2.19E+08	2.19E+08		1	5	0	0	3	1	2	0	3	4	3	2	6
398	ENSG00000136936	XPA	-1	chr9	97687154	97687255		1	64	37	42	48	49	105	27	72	122	111	141	64
399	ENSG00000154767	XPC	-1	chr3	14156348	14156448		1	154	62	86	92	53	241	42	169	230	265	298	178
400	ENSG00000073050	XRCC1	-1	chr19	43552797	43552888		1	40	14	23	21	13	38	8	36	58	50	58	32
401	ENSG00000196419	XRCC6	1	chr22	41621991	41622084		1	940	303	474	450	273	1,188	263	739	1,411	1,435	1,741	886
402	ENSG00000060138	YBX3	-1	chr12	10699501	10699596		1	977	421	429	619	321	1,345	401	1,082	1,523	1,535	1,961	753
403	ENSG00000109906	ZBTB16	1	chr11	1.14E+08	1.14E+08		1	4	2	3	1	3	2	1	0	6	5	6	1
404	ENSG00000128016	ZFP36	1	chr19	39408035	39408161		1	10	3	3	3	2	9	1	5	12	7	5	7
405	ENSG00000185650	ZFP36L1	-1	chr14	68788691	68788787		1	51	20	31	16	22	65	8	32	57	75	81	43
406	ENSG00000162702	ZNF281	-1	chr1	2E+08	2E+08		1	207	99	116	112	67	250	57	190	338	364	415	256
407	_GDC_CONTROL_	_GDC_CO	1	chr3	74811135	74811230	gDNA_control	1	0	0	0	0	0	0	0	0	0	0	0	0
408	_GDC_CONTROL_	_GDC_CO	1	chr1	1.05E+08	1.05E+08	gDNA_control	1	0	0	0	0	0	0	0	0	0	0	0	0
409	_GDC_CONTROL_	_GDC_CO	1	chr21	20087007	20087106	gDNA_control	1	0	1	0	0	0	0	0	0	0	0	1	0
410	_GDC_CONTROL_	_GDC_CO	1	chr18	66284567	66284662	gDNA_control	1	0	0	0	0	0	0	0	0	0	0	0	0
411	_GDC_CONTROL_	_GDC_CO	1	chr15	94718874	94718967	gDNA_control	1	0	0	0	0	0	0	0	0	0	0	0	0
412	_GDC_CONTROL_	_GDC_CO	1	chr15	46219595	46219688	gDNA_control	1	0	0	0	0	0	0	0	0	0	0	0	0
413	ENSG00000130731	C16orf13	-1	chr16	634501	634621	reference_gene	1	113	35	58	50	23	139	36	98	187	158	166	94
414	ENSG00000112787	FBRSL1	1	chr12	1.33E+08	1.33E+08	reference_gene	1	172	76	91	95	52	173	43	143	266	235	326	159
415	ENSG00000100578	KIAA0586	1	chr14	58465870	58465964	reference_gene	1	133	45	68	94	43	146	30	109	162	168	221	114
416	ENSG00000169967	MAP3K2	-1	chr2	1.27E+08	1.27E+08	reference_gene	1	114	58	89	78	46	147	40	121	182	202	270	123
417	ENSG00000084072	PPIE	1	chr1	39743264	39743891	reference_gene	0	199	84	139	140	89	210	66	201	284	284	289	165
418	ENSG00000100023	PPIL2	1	chr22	21684767	21684860	reference_gene	1	205	98	122	107	71	253	57	177	331	324	407	192
419	ENSG00000132005	RFX1	-1	chr19	13961800	13961891	reference_gene	1	54	25	29	21	20	59	10	45	82	74	117	41
420	ENSG00000041988	THAP3	1	chr1	6628529	6628647	reference_gene	1	37	13	14	13	16	36	15	31	44	40	41	25
421	ENSG00000236104	ZBTB22	-1	chr6	33314794	33314894	reference_gene	1	26	20	15	27	15	33	12	32	48	47	37	31
422	ENSG00000083838	ZNF446	1	chr19	58480869	58480974	reference_gene	1	69	20	40	34	12	56	10	53	68	62	75	52

Differential gene expression inter- and intra-samples

Select a Plot or Chart Icon below to launch a separate window.


**Notes:**

1. Before launching a plot or chart, please upload your readout data or take a test run.
2. Please disable any Pop-Up Window Blockers in your web browser.
3. To save the figure:
  - a. Using Windows: Right-click on the figure and save the image/picture
  - b. Using OS X: Hold down the Control key & Click and save the image/picture



Scatter Plot

The scatter plot compares the normalized expression of every gene on the array between two groups by plotting them against one another to quickly visualize large gene expression changes. The central line indicates unchanged gene expression. Set the boundary (fold regulation cut-off) and the experimental groups to compare. Then, export the lists of genes whose expression changes are greater than the selected boundary.



Volcano Plot

The volcano plot helps quickly identify significant gene expression changes. The volcano plot displays statistical significance versus fold-change on the y- and x-axes, respectively. The volcano plot combines a p-value statistical test with the fold regulation change enabling identification of genes with both large and small expression changes that are statistically significant. Select the experimental groups and the boundaries or the fold regulation cutoffs and p-value cutoff values. Then, export the lists of genes whose expression changes and p-values are beyond the selected boundaries.

\* NOTE: This plot requires three or more replicates in each group.



Clustergram

The clustergram performs non-supervised hierarchical clustering of the entire dataset to display a heat map with dendrograms indicating co-regulated genes across groups or individual samples.



Analysis: What kinds of things get flagged?  
 Low tag #, high gDNA, poor normalizer performance




Changes in gene expression due to chemical perturbation and quantified by QIAseq RNA NGS were characterized.

	A	B	C	D	E	F
1	Symbol	Group 1 vs Control Group FC	Group 2 vs Control Group FC	Group 3 vs Control Group FC	Group 4 vs Control Group FC	Group 5 vs Control Group FC
2	ABCC4	0.347303984	0.599495235	0.913120058	1.007402609	0.706931657
3	ABCG2	0.789834324	0.545345791	0.593316022	0.865141275	0.229670718
4	ACACA	0.703476499	0.623639875	0.704593128	0.751523724	0.613677886
5	ACSL3	2.545269175	0.985252008	0.475569118	1.821679828	0.959904487
6	ACSL4	0.607090553	0.954206469	1.16790896	0.920871511	0.783629237
7	ACSL5	1.524483462	0.818680055	8.016242204	1.298761113	3.103061498
8	ACTC1	1.524483462	0.818680055	8.016242204	0.432920371	3.103061498
9	ADA	0.524762261	0.507255299	0.827812466	1.072952252	0.534072559
10	ADCYAP1	1.018793394	0.547113729	5.357155262	0.289315316	2.073737521
11	ADIPOQ	1.524483462	0.818680055	8.016242204	0.432920371	3.103061498
12	ADM	1.287226367	1.814578361	0.846083164	7.265196567	1.637580318
13	ADM2	0.759478555	0.407856143	3.993591403	0.215675503	1.54590635
14	ADORA1	1.524483462	0.818680055	8.016242204	3.896283339	3.103061498
15	AFAP1L2	0.553991157	0.446257558	1.456535084	1.33723293	0.563820032
16	AKR1C2	3.048966924	0.818680055	8.016242204	0.865840742	3.103061498
17	ALAS1	1.862272419	1.072901265	0.998259378	2.259144279	1.490488657
18	ALDOC	5.55235194	1.507248701	0.962509361	7.502548214	0.496779025
19	ANGPTL4	7.097042262	2.177864398	5.331230557	2.015406691	4.127404302
20	ANKRD37	0.712116512	0.600948414	2.139740566	0.866681597	0.828286671
21	APEX1	0.277482896	0.718928251	0.80634339	0.785226601	0.61931136
22	APOA1	0.509396697	0.547113729	2.678577631	0.867945947	1.03686876
23	APOA5	1.524483462	0.818680055	8.016242204	0.432920371	3.103061498
24	APOC3	0.435114402	0.233665691	6.863929684	0.123563091	0.885668346
25	APOE	3.048966924	0.818680055	8.016242204	1.731681484	3.103061498
26	AR	1.018793394	0.547113729	5.357155262	0.578630631	2.073737521
27	ARNT	0.631616725	1.136589464	1.281886459	0.714315259	1.060093787
28	ASNS	2.022391069	1.635229597	0.744168392	1.155112866	4.943570253
29	ATF3	5.781443893	1.395043461	2.464926417	3.610862055	11.4102239
30	ATF4	1.180830418	1.230428098	0.579847466	1.021648438	1.780069371
31	ATM	0.460883341	0.766464495	0.703590455	0.8063953	0.968382561
32	AXIN2	1.161108884	1.402965395	0.381593597	1.154054443	1.329422854
33	BAG2	0.302770298	1.006050942	0.895537771	0.816812128	0.731837318
34	BAX	0.622825729	1.169554298	0.513729458	0.753716486	0.778880146

## Scatter plot and clustergram (HDAC Sample compared to Control)



HDAC is predicted to be inhibited by Trichostatin A and drives a mechanistic network with 18 other regulators.





## QIAseq sample multiplexing guidelines on NGS platforms

			# Samples_ Moderate Read depth			
			average 5000 reads/gene			
			average 2-5 reads per MT			
Instrument	Version	Capacity	100 genes	250 genes	500 genes	1000 genes
MiSeq	V2	15 M	30	12	6	3
MiSeq	V3	25 M	50	20	10	5
NextSeq 500	High Output	400 M	800	320	160	80
NextSeq 500	Mid Output	130 M	260	104	52	26
Ion Torrent PGM	Ion 314 Chip v2	400-550 K	1	0	0	0
Ion Torrent PGM	Ion 316 Chip v2	2-3 M	4	2	1	0
Ion Torrent PGM	Ion 318 Chip v2	4-5.5 M	8	3	2	1
Ion S5	Ion 520 Chip	3-5 M	6	2	1	1
Ion S5	Ion 530 Chip	15-20 M	30	12	6	3
Ion S5	Ion 540 Chip	60-80 M	120	48	24	12

An example: 96 samples, 421 genes

Parameter	QIaseq targeted RNA panels	RT-PCR
Material required	One pool of primers	105 384-well plates
Run time	14 hours for NextSeq run	310 hours (2 hours per plate)
Hands-on time	3 hours (for 96 samples)	105 hours (one hour per plate)
Cost per sample	\$65 (exclusive of sequencing run)	\$239
Sample	10 ng each sample	4 µg each sample






Molecular barcoding allows for traceability of variants, AND CNV analysis.




QIASeq Targeted RNAscan is a RNA target enrichment method that allows verification of known fusions and discovery of novel fusions with next-generation sequencing (NGS).

1. [NCOA3--SULF2](#), score: 3806.00
2. [RPS6KB1--VMP1](#), score: 1136.00
3. [ARFGEF2--SULF2](#), score: 184.00
4. [SULF2--PRICKLE2](#), score: 56.00
5. [CLTC--VMP1](#), score: 6.00

## RPS6KB1-VMP1




## ARFGEF2-SULF2





- 1 Using NGS approaches for gene expression analysis
- 2 Principle of QIAseq Targeted RNAseq
  - 2.1 Molecular Barcodes
  - 2.2 QIAseq RNA workflow
- 3 An application of the QIAseq RNA system
  - 3.1 QIAseq data analysis
  - 3.2 Ingenuity IPA
- 4 Uses of UMIs in other types of targeted sequencing**
- 5 Summary and Discussion

# QIAseq miRNA Sequencing Kit: One-day Workflow




● Sample to sequencer in 7 hours! Other protocols can take **2 days!**



- 1 Using NGS approaches for gene expression analysis
- 2 Principle of QIAseq Targeted RNAseq
  - 2.1 Molecular Barcodes
  - 2.2 QIAseq RNA workflow
  - 2.3 QIAseq RNA performance
- 3 An application of the QIAseq RNA system
  - 3.1 QIAseq data analysis
  - 3.2 Ingenuity IPA
- 4 Summary and Discussion**

- Extremely sensitive expression profiling, >1 copy per cell
- Highly flexible in experimental design, from 12 to 1000 or more targets, 1 to 96 samples
- High specificity, ~97-99% maintained through all panels
- Extremely high primer uniformity ~0.98 at 20% mean
- Random molecular barcoding for quantification
- Requires no rRNA depletion or blocking or dT selection
  - Only requires ~1ng-20ng total RNA
- Makes best use of limited NGS read budget
- System optimized for best possible performance with FFPE samples
- Leverage QIAGEN content know-how for NGS
  - Disease and pathway specific collections
  - Extended panels and fully custom gene content 12-1000 genes





# QIAseq Targeted RNA Products

## QIAseq Targeted RNA Panel (12 or 96 samples)

Kit containing reagents for first strand synthesis, Smcounter tagging, and gene-specific amplification for targeted RNA sequencing

## QIAseq Targeted RNA Extended Panel (12 or 96 samples) (up to 25 additional targets)

Kit containing reagents for first strand synthesis, Smcounter tagging, and gene-specific amplification for targeted RNA sequencing;

## QIAseq Targeted RNA Custom Panel (12, 96 or 384 samples)

Kit containing reagents for first strand synthesis, Smcounter tagging, and gene-specific amplification for targeted RNA sequencing

## QIAseq Targeted RNA sample Indexing(12-plex or 96-plex HT) for ion torrent

## QIAseq Targeted RNA sample Indexing (12-plex or 96-plex or HT) for illumina

## Library Quant Assay/Array Kit

Assays and master mix for library quantification prior to NGS

---

### **Initial Content:comprehensive 250 – 500 gene panels and ALL human RT2 panel content (200 panels)**

Immunity and Inflammation

Angiogenesis and Endothelial

Cell Death

Cancer Pathway

Signal Transduction

ECM and Cell Adhesion

Molecular Toxicology

Stem Cells

# Questions?



## Contact QIAGEN Technical Service

Call: 1-800-426-8157 for US

Call: +49 2103-29-12400 for EU

### Email:

[DigitalRNAseq@qiagen.com](mailto:DigitalRNAseq@qiagen.com)

[techservice-na@qiagen.com](mailto:techservice-na@qiagen.com)

[techservice-eu@qiagen.com](mailto:techservice-eu@qiagen.com)

Webinar related:

[QIAwebinars@qiagen.com](mailto:QIAwebinars@qiagen.com)