



Human ID and Forensics

Get the most out of every sample

With ultra-efficient nucleic acid extraction on EZ2[®] Connect Fx



EZ2 Connect Fx Instrument

Simplicity. Safety. Speed.

Designed with forensics labs in mind

Fast and efficient workflow


- Purify high-quality DNA and RNA in under 20 minutes
- Increase lab throughput by processing up to 24 samples in parallel and 288 samples per day – just with one instrument!

Maximum process safety

- Achieve reproducibility and convenience with prefilled, ready-to-use forensic grade (ISO 18385) reagent cartridges
- Have confidence in your sample processing with an internal camera for load checks and sample recovery capabilities
- Full worktable decontamination in one easy step with high emission UV-LED

Trusted chemistry

- Benefit from the EZ1® Advanced XL magnetic-bead purification technology you know and trust
- Easily access a variety of dedicated human ID and forensics protocols at the push of a button
- Prepare a wide range of sample types for various downstream forensics applications



Many sample types



Trace samples

Body fluids





Bone

Tissue



Standard samples

Other human ID samples



Many applications

PCR




qPCR



NGS



Get high-quality nucleic acids in just 20 minutes



Visit www.qiagen.com/HID-EZ2-connect-Fx

and unravel the simplicity, process safety and speed in HID sample prep

Automate your workflow

Confidence and convenience are now within your reach.

Contact your sales representative

to find out how you can get an EZ2 Connect Fx for your lab.



The EZ2 Connect Fx is intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

Trademarks: QIAGEN®, Sample to Insight®, QIASphere®, EZ1®, EZ2® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. PROM-19447-001-1126213 10/2021 © 2021 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com